

## VM-80 Vertical Rotating Mixer

VM-80 Vertical Rotating Mixer provides efficient but soft mixing to keep the sample in a suspension status such as blood mixing. It is suitable for preventing blood clots, latex diagnosis, immunoprecipitation and similar applications.



### 【 Product Features 】

1. Smart appearance, solid structure, simple operation. LED display time and speed.
2. DC brushless motor, stable performance, low noise and long using life.
3. 360° vertical circulation rotation, operating smooth and steady.
4. Different tube holder frames to option to suit for different kinds tubes.

### 【 Technical Parameters 】

Model	VM-80
Speed Range	10rpm ~80rpm
Shaking Orbit	circle
Time Range	1min ~99h59min or continuous
Voltage	AC 100 ~230V, 50/60Hz
Power	16W
Fuse	250V, 1A, $\Phi$ 5x20
Dimension (WxDxH)	470 x 190 x 160mm
Net Weight	2.8kgs



360 Rotation Angle



LED Display



DC Brushless Motor



3 Platforms (platform x 3)



### 【 Optional Accessories 】



VM80-A Platform: 1.5ml Centrifuge x 64



VM80-B Platform: 15ml Sharp Bottom Tube x 22



VM80-C Platform: 50ml Sharp Bottom Tube x 16