



OS-20 ILLUSTRATED WITH
PP-4 PLATFORM

BOECO UNIVERSAL ORBITAL SHAKER OS-20

Compact, low-profile, orbital laboratory shaker for versatile applications in the fields of biotechnology, microbiology and medical diagnostics. The OS-20 provides gentle rotational mixing and precise circulation at each point of the platform. Variable digital displayed speed control from 50-250 rpm. Digital timer for setting a time interval and continuous operation.

Specification:


Kind of motion:	orbital
Orbit size:	20 mm
Speed regulation:	50 - 250 rpm
Speed display:	digital (LED)
Timer:	from 1 min to 999 min / continuous
Time display:	digital (LED)
Max. load:	2,5 kg excl. attachment
Power:	12 V. DC 2,0A (External power adapter)
Overall size :	265 x 270 x 70 mm (w/o attachment)
Weight:	3,2 kg.

Code

BOE 8059000

Description

Universal Orbital Shaker OS-20, AC Adapter 100-240V with Multi plug (EU, US, UK, AU)



PSU-10i ILLUSTRATED WITH
UP-12 PLATFORM

BOECO UNIVERSAL ORBITAL SHAKER PSU-10i

Equipped with the direct drive system and brushless motor, our compact, powerful and quiet microprocessor controlled orbital laboratory shaker for many general purpose shaking applications. Speed and time setting are easy and during operation the LCD status display clearly indicates both set and actual speed, plus set and elapsed time.

Specification:

Kind of motion:	orbital
Orbit size:	10 mm
Speed regulation:	50 - 450 rpm
Speed display:	digital (LCD)
Timer:	from 1 min to 96 hours / continuous
Time display:	digital (LCD)
Max. load:	3 kg excl. attachment
Power:	12 V. DC 2,0A (External power adapter)
Overall size:	220 x 205 x 90 mm (w/o attachment)
Weight:	3,4 kg

Code

BOE 8069200

Description

Universal Orbital Shaker PSU-10i, AC Adapter 100-240V with Multi plug (EU, US, UK, AU)

ACCESSORIES FOR OS-20 AND PSU-10i

PLATFORM PP-4

Code

BOE 8059001

Description

Dish platform PP-4 with silicone mat for Petri dishes, culture flasks etc. 235 x 235 mm

PLATFORM UP-12

BOE 8059002

Universal platform UP-12 for all kind of flasks 285 x 215 mm, stainless steel with variable clamping rolls and silicone mat.

PLATFORM P-6/250

BOE 8059003

Numbers of flasks at the platform
100- 150 ml flasks max 12-9 pieces
200- 250 ml flasks max 6-4 pieces
500-1000 ml flasks max 4-2 pieces

PLATFORM P-12/100

BOE 8059004

Platform P-6/250 with fixed clamps for 6 flasks 250/300 ml, 250 x 190 mm
Platform P-12/100 with fixed clamps for 12 flasks 100/150 ml, 250 x 190 mm

