

VM-80 Vertical Rotating Mixer


Rotator, Mixer and Vortexer

Introduction

Vertical Rotating Mixer provides efficient but soft mixing to keep the sample in a suspension status such as blood mixing. It is suitable for preventing blood clots, latex diagnosis, immunoprecipitation and similar applications.

Feature

1. Smart appearance, solid structure, simple operation. LED displays time and speed.
2. DC brushless motor, stable performance, low noise and long using life.
3. 360°vertical circulation rotation, operating smoothly and steadily.
4. Different tube holder frames for optional to suit for different kinds tubes.
5. Oil-free bearings are built in both ends of the bracket.



Parameter

Model	VM-80
Speed Range	10rpm~80rpm
Shaking Orbit	circle
Time Range	1min~99h59min
Voltage	AC 100~230V, 50/60Hz
Power	16W
Fuse	250V, 1A, Φ 5x20
Dimension	W.470 x D.190 x H.160mm
Net Weight	2.8kgs

Accessories



VM80-A Tube Holder Frame:
1.5ml Centrifuge x 64



VM80-B Tube Holder Frame:
15ml Sharp Bottom Tube x 22



VM80-C Tube Holder Frame:
50ml Sharp Bottom Tube x 16